

Technical Datasheet

		EN ISO 13485 certified
Article - No:	07-902-0300	Revision: 2
Remarks		single use only

Specification

Product	Gel-loading pipette tip PP, universal fit, 0,1-10µl, round tip, 200 pcs/rack, transparent, np pcr ready
Dimensions	length: 45,90 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	0,1-10µl
Material Characteristics	polypropylene (PP)
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per inner cs: 200
	pcs per carton: 800
	pcs per pallet: 288.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	